

# MAGNETIC PULL CORD

The magnetic pull cord is great for common areas, restrooms, as well as resident rooms and bathrooms.



PAL -203200

## Features

- Perfect for common area bathrooms, resident rooms, bathrooms, and common areas.
- Wireless transmitter communicates with PalCare server.
- Long Life 3V Lithium Battery (replaceable on site).
- Call is activated when pulled, and can be reset by lifting up on the emergency switch.
- Available in surface mount or flush mount.

## Specifications

- Dimensions: 5" high, 3.25" wide, 2.25" deep (surface mount box).
- Color: White with Red Emergency Call Switch.
- Power: CR123A lithium battery
- Battery Life: 1 year @ 3 activations / day
- CleanCord™ for easy cleaning and sterilization
- Compatibility: Works with all PalCare nurse call systems.